geat2GO indigo. Tailored solutions. Greater independence.

Product Information Guide

Guardian Mediwatch Community Based/Mobile Watch Alarm

(Supplier: Guardian)



| Monitoring | Both family monitoring and professional monitoring options available | | |
|--|---|--|--|
| Ongoing Costs (correct at date below, subject to change as per suppliers' advice) | Family monitored: Optus \$50/year; Telstra \$70/year (default); TPlus (Unlimited calls) \$150/year (<i>first year covered by geat2GO</i>) Monitored: \$31/month (first 2 months covered) | | |
| Colour | Black | | |
| Case material | Plastic | | |
| Accessories | Magnetic charger and silicone wrist band (included), Other wrist band options available via Guardian directly) | | |
| Range/Reception | Community use with mobile coverage (2-3 bars reception required) | | |
| Water resistance | Shower/splash proof (avoid submerging) | | |
| Dimensions (H x W x D) | 44 x 46 x 17mm | | |
| Weight | 63g | | |
| Battery life | Up to 3 days depending on usage (daily charging recommended) | | |
| Fall detection | Yes | | |
| Live tracking/GPS | Yes | | |
| Two-way call | Yes | | |
| Geofence | No | | |
| No Movement Detector Alarm | Yes | | |
| Heart rate monitor | Yes (within 10 bpm) | | |

indigosolutions.org.au

© 2023 Indigo Australasia Incorporated

ice.

| geat2GO | | inc | lig | 0. | Tailored solutions. Greater independenc |
|---------|----------|-----|-----|----|--|
| | . | | | 14 | |

| Step counter | Disabled by default, enabled by request | | |
|------------------------------------|---|--|--|
| Voice Prompts, Languages | Yes, English included or other languages available at \$30 additional cost | | |
| Pacemaker information | Consult your cardiologist/specialist prior to using the device if you have a pacemaker. | | |
| Emergency contacts | Can have up to 10 contacts for family monitoring | | |
| Set up needed? | Connect charger to power | | |
| SOS procedure | Press and hold the SOS button until it vibrates. A voice prompt will then announce the activation. Family Monitored: A call will be made to Contact Number. All emergency contacts will receive an SMS. Professionally Monitored: The Monitoring Company will receive an alert and call the MediWatch within a few moments. After the call, the MediWatch will announce a reminder to press the SOS button which will end the SOS sequence, meaning that no more emergency contacts will be called. | | |
| How to cancel alarm | Press the X on the screen within 10 secs, or at any time | | |
| Additional settings/information | Settings menu disabled by default, can request to be activated (Settings: Power off/restart, change watch face, view watch information) Watch-face able to be changed between digital and analog within watch settings. SIM whitelisted by default | | |

Information correct as at 30 October 2023